


ISPE-CASA PRESENTS “CLEAN & GREEN”



CLEAN AND GREEN TECHNOLOGIES FOR PURE WATER AND PROCESS TECHNOLOGIES




Sequential Multi-Column Chromatography

ISPE, April 13, 2009

Michael LaBreck
Sales Manager
Novasep, Inc.




Presentation Outline

- Purpose and Applications of SMCC
 - Functional Principles
 - Case Study including Economic Evaluation
 - Equipment Product Specification
 - Practical & Regulatory Considerations
 - Conclusion
- 



Sequential Multi-Column Chromatography : Purpose & Application

- “Debottlenecking” Downstream Purification Processes
 - **Improving Efficiency of all Chromatography Applications**
 - Capture or Flow Through
 - Affinity, IEX, HIC; practically all types of stationary phases
 - Applicable to Low, Medium and High Pressure Chromatography
 - **Applicable to all Types of Biomolecules**
 - MAb's
 - rTherapeutics
 - Peptides
 - Blood Fractions
 - Vaccines
- 

Current Chromatography Productivity Paradigm


Affinity Capture Example

Batch
Affinity Chromatography




Increased
Loading Capacity
(50 mg/ml – 100 cm/h)

OR



Increased Feed
Flow
(24 mg/ml – 450 cm/h)

SMCC
Affinity Chromatography



Increased Loading
Capacity


AND

Increased Feed Flow
(50 mg/ml – 450 cm/h)

- Thermodynamics “Load/Resin interactions”
Vs Kinetics “Flow Rate”
- Minimum process flow rate requirements compromise maximization of productivity
- Significantly higher loading flow rates
- Loading on multiple smaller columns
- Minimizes impact of kinetics ⇒
Increased dynamic resin capacity (g/L)

Sequential Multi-Column Chromatography

Functional Principles



Sequential Multi-Column Chromatography Batch / SMCC Comparison

Batch Chromatography

- **Wasted Resin Capacity**
 - Resin capacity not used 100%
- **Wasted Water**
 - BV's for wash, elute, regeneration, etc based on actual bed volume, not working bed volume
- **Wasted Time**
 - Top layer is waiting to be eluted
 - Lower layer is waiting to be loaded

SMCC


- **Less Resin**
 - 3-6 Smaller Ø columns
 - Typically >40% resin volume reduction
 - Smaller space requirements
- **Less Water**
 - Smaller bed volumes
 - Same BV elute ratio
 - Less buffers
- **Faster Processing Time**
 - Higher productivity

Sequential Multi-Column Chromatography


Cycle on 1 column



Chronogram




Batch Process Chromatogram




Sequential Multi-Column Chromatography

Batch Characterization 100 to 600cm/hr

Breakthrough Curves



Loading Capacity




Sequential Multi-Column Chromatography

Loading Capacity

Multicolumn Process



Single Batch Column @ 100cm/hr



Sequential Multi-Column Chromatography

Loading Capacity

Multicolumn Process



Single Batch Column @ 100cm/hr




Sequential Multi-Column Chromatography

Loading Capacity

Multicolumn Process



Single Batch Column @ 100cm/hr



BioSC -2 Column SMCC Process @ 100cm/hr

Sequential Multi-Column Chromatography

Loading Capacity

Multicolumn Process



Single Batch Column @ 100cm/hr




BioSC -2 Column SMCC Process @ 100cm/hr

Sequential Multi-Column Chromatography

Column Loading Capacity Profiles



Single Batch Column @ 200cm/hr





BioSC -2 Column SMCC Process @ 200cm/hr

Sequential Multi-Column Chromatography

Column Loading Capacity Profiles



Single Batch Column @ 400cm/hr




BioSC -2 Column SMCC Process @ 400cm/hr

BioSC -3 Column SMCC Process @ 400cm/hr



Sequential Multi-Column Chromatography

BioSC 4 Column SMCC Flow Schematic - 1/4




Sequential Multi-Column Chromatography

BioSC 4 Column SMCC Flow Schematic - 2/4




Sequential Multi-Column Chromatography

BioSC 4 Column SMCC Flow Schematic - 3/4



Sequential Multi-Column Chromatography

BioSC 4 Column SMCC Flow Schematic - 4/4



Sequential Multi-Column Chromatography

Lab Scale Affinity Capture – Mab @ 1g/L

		Experimental Data		Simulation
		Single Column Batch	BioSC 3 Column	BioSC Simulation
# Columns		1	3	3
Resin Volume (ml)		60	60	60
Dynamic Capacity (g mAb/L resin)	300cm/hr	25 g/L		40 g/L
Buffer Consumption (ml/g mAb)		1200 ml/g		675 ml/g
Productivity (g mAb/L resin/day)		200 g/l/day		215 g/l/day
Dynamic Capacity (g mAb/L resin)	450cm/hr	20 g/L	40 g/L	40 g/L
Buffer Consumption (ml/g mAb)		1500 ml/g	830 ml/g	675 ml/g
Productivity (g mAb/L resin/day)		220 g/l/day	310 g/l/day	320 g/l/day
Dynamic Capacity (g mAb/L resin)	600cm/hr	15 g/L	40 g/L	40 g/L
Buffer Consumption (ml/g mAb)		2000 ml/g	675 ml/g	685 ml/g
Productivity (g mAb/L resin/day)		240 g/l/day	430 g/l/day	445 g/l/day

- o 44% or 1.8 fold increase in productivity
- o 2.9 Fold decrease in water usage
- o Purity 99% in each case
- o Yield 98% in each case

- o Excellent agreement between simulation modeling & empirical data
- o Implications for development & optimization of continuous processes

Commercial Protein A Capture Process – Simulation @ Const. Resin Volume

		Batch mAb capture	BioSC SMCC Modelling Data		BioSC SMCC Modelling Data
		Single Column Batch Process	3 Columns	4 Columns	4 Columns
Number of columns		1	3	4	4
Bed Height (cm)		20	20	20	14
300 cm/h	Dynamic capacity (g mAb/L resin)	16.8 g/L			
	Buffer Consumption (mL/g mAb)	1500 mL/g			
	Productivity (g mAb/L resin/day)	386 g/L/day			
750 cm/h	Q'max Capacity (g mAb/L resin)		40.8 g/L		
	Buffer Consumption (mL/g mAb)		630 mL/g		
	Productivity (g mAb/L resin/day)		786 g/L/day		
1200 cm/h	Q'max Capacity (g mAb/L resin)			40.8 g/L	
	Buffer Consumption (mL/g mAb)			630 mL/g	
	Productivity (g mAb/L resin/day)			939 g/L/day	
950 cm/h	Q'max Capacity (g mAb/L resin)				40.8 g/L
	Buffer Consumption (mL/g mAb)				630 mL/g
	Productivity (g mAb/L resin/day)				969 g/L/day

• **Best Scenario:**

- 4 Columns with resin volume as batch process
- Smaller diameter columns
- 2.6 Fold increase in productivity
- 2.4 Fold reduction in buffer usage

40.8 g/L
630 mL/g
969 g/L/day

Generic Commercial Protein A Capture Step

Best Batch Scenario Single Step Purify 10 kg-1 Cycle

L/Resin/day for kg at 386 g/L	26 L
L/Resin for 10 kg mAbs	624 L
Column Size for 20 cm BH	2.0 meter
Daily production	10 kg
Assuming 200 Regeneration Cycles	3.4 kg per lifetime Liter

SMCC Scenario Single Step Purification 10 kg-1 Cycle

L/Resin/day for kg at 969 g/L	10 L
L/Resin for 10 kg mAbs	250 L
Column Size (4 columns with 14 cm BH)	75 cm
Daily production	10 kg
Assuming 200 Regeneration Cycles	8.0 kg per lifetime Liter

Economic Implications of SMCC process

Batch - 1.7meter diameter with 20cm bed height, DBC 20mg/ml, 300 cm/hr

SMMC – (4) 80cm diameter with 10cm bed height, DBC 45 mg/ml, 600cm/hr

	Batch MAb Capture	SMCC MAb Capture	Potential Savings
Volume of Protein A Media	455L	200L	>\$3,000K
Process Time	4 hours	2 hours	50% savings
Buffer Usage	18,000 L	5,800 L	66% savings

BioSC™ Hardware Features

BioSC™ Pilot System

Capacity & Equipment

- 6 x 5 - 15 cm columns (Glass, Acrylic or Stainless Steel)
- In-Line Buffer Dilution
- Isocratic, Gradient or Step gradients
- Analog «mass flow meter» controlled
- Up to 4 fraction collections as standard
- Up to 8 inlet, load, buffer lines
- Fully Automated
- SIP option
- On Line/At Line Monitoring
- Footprint : 2 m x 2 m x 1 m

ASME/BPE Engineering Design

Materials:

- 316 L Stainless Steel
- Teflon, EPDM, Glass or Acrylic
- Sanitary Design for Assured CIP
- Diaphragm valve
- Rotating piston pumps
- >1% pipe slopes
- Fully drainable
- No dead legs
- Column bypass loops
- TC connections

Control Platform


- Chromsoft user friendly interface
- Intellutions iFix control systems
- 21 CFR Part 11 compliant
- GAMP 4 compliant



- *Small Production Foot Print*
 - 2m x 2m x 1m
- *LabScale System Available*
 - 0.5 – 2.6 cm Ø Columns
- *Process Scale*
 - Specification to Client URS



Regulatory and Quality Considerations for SMCC

- **Regulatory and Quality Questions:**
 - Still Chromatography using the same buffers ...
 - **What is a batch?**
 - Typically defined by the process strategy, by the starting material (e.g. the fermentation batch) or by the formulation quantity and even sub lots.
 - The Overall DSP process remains batch even if 1 chromatography step is changed to sequential.
 - **On Line Control:**
 - Critical quality attributes can be controlled On or At Line with BioSC™, compliant with PAT initiatives.
 - **Quality Profile of a continuous or “steady state” process is much more homogenous than batch**
 - Stationary phase usage is more homogenous.
 - Less overloading, less fouling, less cleaning.
 - Improved performance of chromatography resin.
 - Scalable process
- 



Regulatory and Quality Considerations for SMCC

- **Typical Strategy for SMCC Process Validation**
 - Characterization of the batch process (Breakthrough curves with the purified molecule & the load)
 - Process Simulation
 - Lab tests to verify the predicted process
 - Resin life study
 - Test of the robustness / definition of PAR of the process -> Design Space
- **Implementation at the industrial scale including IQ, OQ, PQ of equipment**
- **Process Validation**

Note: SMCC uses the same buffers, CIP solutions and stationary phases as batch processes



Regulatory and Quality Considerations for SMCC


FDA Approved API Process

- Continuous, fully automated systems
- DAC columns
- 1 meter in diameter
- Approved by FDA
 - **Varicol**®






SMCC Processing Using BioSC™ In Conclusion

- **Established scientific principles**
 - Applicable to most chromatography types
 - **Utilizes same product contact materials**
 - Resins, Buffers, Eluents, etc.
 - **Delivers significant economy of scale**
 - 2-4 Fold increase in productivity
 - 2-4 Fold reduction in water usage
 - Positively impacts facility CapEx, foot print, utilities, tanks, etc.
 - 2-4 Fold increase in process speed
 - **Manageable implementation and learning curve**
 - NOVASEP Expertise
 - Running FDA regulated continuous processes
 - Simulation & modeling capability
 - NOVASEP Biopharma PPD™ expertise for contract process development
- 



Acknowledgements

- **Margit Holzer, Laurent David, Jean Guillerm & H. Osuna-Sanchez**
Novasep Process
- 



Thank you for your attention!

**Michael LaBreck
Novasep, Inc.**

michael.labreck@novasep.com